



FDA EUA AUTHORIZATION  
OVER THE COUNTER

CareStart™

COVID-19

ANTIGEN  
HOME TEST

## CareStart™ | ANTIGEN HOME TEST

2 Test Per User Box

### Features:

1. FDA EUA for at HOME USE
2. Lateral flow immunochromatographic assay
3. Using anterior nasal swab (appropriate for self-collection)
4. Test results are interpreted at 10 minutes after sample loading
5. Visual Read or App Read for result interpretation (no instrument required)
6. Digital guide for self sample collection and self-testing (digital IFU in mobile application)
7. 2D-barcode on each device: ensuring traceability of device (manufacturing information)
8. Digital connectivity: automatic result reporting to local or federal health authorities

### Palletized // Weights & Measures

6400 - 2 Pack Boxes Per Pallet

48.00 x 40.00 x 75.50

170 KG

**TRACKS ALL KNOWN VARIANTS**

*Per Manufacturer's June 10 Letter Available Upon Request*